


∧∨∣PURE®-WT

HIGHLY RESISTANT PTFE INLINER WITH SILICONE COVER AND STAINLESS STEEL WIRE



PRODUCT INFORMATION

AVIPURE® - WT is designed to facilitate pressurized fluid transfer in the biopharmaceutical manufacturing processes

AVIPURE® - WT has an inliner made of Teflon® PTFE, which ensures a highly chemical resistance. The cover design made of silicone with braid and stainless steel wire reinforcement gives a high pressure resistance and maintain an outstanding flexibility.

The AVIPURE® - WT is available in 6 different dimensions. More sizes on request.

APPLICATIONS

AVIPURE® - WT is designed to be used in many pharmaceutical and biopharmaceutical applications where chemical, pressure and vacuum resistance and a good bending radius is required such as:

- → Transport of liquid or semi-liquid fluids
- → Vessel or tank transfer
- → CIP solution

FEATURES AND BENEFITS

- → "Extractable studies" for AVIPURE® WT available
- Pure PTFE: Resistant against acids and bases and extraordinary heat and cold resistance
- Continuous extrusion of inliner allows smooth internal surface and optimum cleaning results
- → Small bending radius
- Marking of product information on the hose for perfect treacablity
- → Coils wrapped in double PE-bags
- "Low-Tack" significantly reduced surface friction:
 Easier handling with gloves
- → Certificate of Conformance: Delivered with each order for traceability

FLEXIBILITY

AVIPURE® - WT is available as hose coils of 10 m. Cutting of hose without special tools. Fits perfectly to the reusable hose connection AVILOC®.



Technical contents are subject to change

∧∨∣PURE®-WT

HIGHLY RESISTANT PTFE INLINER WITH SILICONE COVER AND STAINLESS STEEL WIRE



TECHNICAL DATA

SPECIFICATIONS:

Colour: Inliner white, cover: natural silicone Material: Teflon® Inliner with Silicone cover Shore A hardness: 60 ± 5 (Silicone)

Temperature range: from -40°C to +150°C

Tear strength: > 9.0 MPa **Elongation at break:** > 400% Compliant: USP <88> Class VI, FDA

Marking: Layline with Type, Size, Lot-No and MFD

Sterilization: Autoclave, SIP

Hose coils packaging: Single PE bag

Shelf life: 5 years

AVIPURE® - WT IS TESTED FOR COMPLIANCE TO:

- **Teflon® PTFE Inliner:**
- FDA 21 CFR 177.1550
- USP <88>: Class VI Implant test, systemic toxicity tests and intracutaneous tests
- Cytotoxicity according to ISO 10993-5
- EU1935/2004 EU10/2011
- 3A Sanitary Standard
- REACH

TYPE	ID X OD DIMENSI- ONS	VACUUM	WORKING PRESSURE AT 20°C		MINIMUM BEND RADIUS	HOSE COIL LENGTH	WEIGHT	ARTICLE- NUMBER
	MM	BAR	BAR	BAR	MM	М	KG/M	
WT 500	13 × 24	0.9	10	40	60	10	0.4	407HSWT0500-10
WT 750	19 × 30	0.9	10	40	90	10	0.7	407HSWT0750-10
WT 1000	25 × 36	0.9	10	40	115	10	0.8	407HSWT1000-10
WT 1250	32 × 43	0.9	8	32	150	10	1.1	407HSWT1250-10
WT 1500	38 × 50	0.9	7	28	180	10	1.4	407HSWT1500-10
WT 2000	50 × 62	0.9	7	28	230	10	1.8	407HSWT2000-10

Other sizes on request

TYPE	ID X OD DIMENSI- ONS	VACUUM	WORKING PRESSURE AT 68°F		MINIMUM BEND RADIUS	HOSE COIL LENGTH	WEIGHT	ARTICLE- NUMBER
	INCH	PSI	PSI	PSI	INCH	FT	LBS/FT	
WT 500	0.5 x 0.94	13	150	580	2.36	32.8	0.28	407HSWT0500-10
WT 750	0.75 x 1.18	13	150	580	3.54	32.8	0.46	407HSWT0750-10
WT 1000	1.0 x 1.42	13	150	580	4.53	32.8	0.54	407HSWT1000-10
WT 1250	1.25 x 1.69	13	115	460	5.91	32.8	0.74	407HSWT1250-10
WT 1500	1.5 x 1.97	13	105	410	7.09	32.8	0.94	407HSWT1500-10
WT 2000	1.97 x 2.44	13	105	410	9.06	32.8	1.21	407HSWT2000-10

Other sizes on request

- The bursting pressure of the AVIPURE* WT hose is reduced by around 20% per temperature increase of 93°C until the maximum hose temperature is reached. The reduction is linear.
- Pressure limit rates can be limited by used end connections.
- · Norm: ISO 1307 for dimensional tolerances

Technical contents are subject to change

